















CUSTOMER PRODUCT SPECIFICATION

Page 1 / 2

Product code	8.01	Version	1.2	Holder	DB
Drawing nr.	8.01_v1.0	Date	6-10-2020	Approval	

Product code	8.01	Sterilization process	EO (Ethylene Oxide)
Description	Self-adhesive ETO dot.		
Validated range	Parameters: > 25 min >37° C >600 mg/l 100% ETO// or >20 min, >54° C >600 mg/l. The colour change is not limited to afore mentioned parameters, but actual colour change has to be determined		
Colour change	Violet to green.		
Application	This process indicator is designed for use with individual items (e.g. packs, containers) to show that the unit has been directly exposed to the sterilization process and to distinguish processed from unprocessed items.		


Size (w x h)	Ø 12,7 mm	Roll diameter	Max 125 mm
Label material	Self-adhesive wood free paper	Shelf life	24 months after date of production.
Adhesive	Permanent	Storage	Dark and dry, between 10° and 30°C or 50° and 86°F.
Roll width	14,2 mm	Imprints	-- -- --
Quantity/Roll	5000 labels	Colours	-- -- --
Core: Ø / material	46,0 mm Cardboard	PMS-nrs.	-- -- --

Winding nr.	-								
	-								
	Inside	 Outside				 Inside			

Packaging	10 rolls packed in a cardboard box.
Product label	Lot # (label) / Productcode / Minimum dose / Colour change / Shelf life / Quantity / Storage.
Remarks	--

Product code	8.01	Version	1.2	Holder	DB
Drawing nr.	8.01_v1.0	Date	6-10-2020	Approval	

Drawing	8.01_v1.0
Date	6-12-2012
Scale	Not drawn to scale.
Colours	Screen colours can differ from actual product colours.


Customer approval

Customer:

Signature :

Date :

Name :