

**WarwickSASCo**  
Moulding The Future of Medical Plastics



## 5 reasons to choose Warwick SASCo

Recognised as the market leader in the design and  
manufacture of medical polyware

Wide range of high quality re-usable and single-use polypropylene  
products for all areas of healthcare

Unrivalled reputation for providing quality products, delivery and  
service to healthcare professionals worldwide

Re-usable products are guaranteed for a minimum of 1000 wash  
cycles – delivering excellent value for money

Certification since 1997 demonstrates our continued commitment  
to maintaining industry standards. Our quality management system  
complies to ISO1345:2016. Products comply with 93/42/EEC  
Medical Devices Directive Class 1 and show CE where applicable.

 **MADE IN BRITAIN**



## Warwick SASCo

For 35 years Warwick SASCo has been part of evolving supply of  
decontamination and nursing products worldwide to hospitals, healthcare and  
patient nursing.

Recognised as a market leader, we have a reputation for high quality products  
designed for hospital decontamination and sterile services, operating and  
surgical applications, patient care and home nursing - where infection  
management and good user application is paramount.

Warwick SASCo's low cost reusable and single use product ranges are chosen  
by healthcare professionals and end users worldwide.

Our growing family owned business is proud to be represented by leading  
distributor partners in over 70 countries.

All our products are manufactured within the UK with Quality Assured  
Suppliers. Our customer support and distribution centre is in Warwick, England.

The team at Warwick SASCo is immensely proud of delivering high quality  
products and a professional friendly service to you - our customer. You will find  
our experienced team are knowledgeable, caring and quick to respond to  
your requirements.

Darby Booth,  
Managing Director – Warwick SASCo



# Products Overview

Our wide range of reusable and disposable medical polyware products are chosen by healthcare professionals for their high quality and low cost.

## Reusable Medical Polyware

WARWICK RANGE

The Warwick Range for use in Decontamination Services, Operating Rooms, Surgeries, Clinics and healthcare services. Products are designed for repeated washing, disinfection and sterilization.

The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>). Products comply with 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.

## Single Use Medical Polyware

STERILITE RANGE

Sterilizable lightweight polyware for single-use and disposable sets. SteriLite Range is designed for Custom Pack Manufacturers and Decontamination Services.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>). Products comply with 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.

## Nursing & Care Daily Living Aids & Accessories

SPA RANGE

Spa Range of personal living aids, patient care, toileting and bedside products.

Products comply with BS2588 For Portable Sanitary Products, BS3215 for Reusable Urinals, 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.



## How to Order

To place an order, request a quotation or for general enquiries:

CALL us on +44 (0)1926 422 427  
FAX +44 (0)1926 315 405  
EMAIL [enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

We work with a wide range of industry sectors including:



### STERILE SERVICES & DECONTAMINATION

Operating Room, Instrument Packing



### PERSONAL USE & HOME CARE

Daily Living Aids, Community Nurse, Incontinent Care



### HOSPITALS

Surgery, Clinics, Wards, Veterinary, Dentistry



### NURSING HOMES

Patient Care, Rehabilitation



### CUSTOM PACK MANUFACTURERS

Single Use Packs, Contract Instrument Processing

## WARWICK - reusable range

Compartment Trays	page	6
DIN Style Baskets	page	7-8
Disinfectant & Soaking Containers	page	6
Gallipots	page	9-10
Instrument Trays	page	7-8
Instrument Tray Lids	page	9
Jugs	page	10
Kidney Dishes	page	11
Lotion & Sponge Bowls	page	9
Medical Boxes, Instrument Boxes & Storage Containers	page	12
Scope Containers	page	13
Silicone Protection	page	13
Theatre Bowls	page	11
Tray Tags	page	14
Wash Bowls	page	11

## STERILITE - single use range

Compartment Trays	page	16
Gallipots	page	17
Jugs	page	18
Kidney Dishes	page	18
Packing Trays	page	16
Quivers	page	18
Sponge Bowls	page	17
Theatre Bowls	page	18

## SPA - care and living aid range

Bedpans	page	20-21
Bedpan Holders	page	21
Bidets	page	21
Commode Pans	page	21
Denture Cups	page	22
Drinking Beakers and Cups	page	22
Drinking Cup Lids	page	22
Jugs	page	24
Jug Sets	page	24
Kidney Dishes	page	23
Medicine Measures	page	24
Patient Trays	page	23
Urinal Bottles	page	25
Urinal Pans	page	25
Vomit Bowls	page	23



To find out more about our service, product ranges, product decontamination cycle and industry news:

visit [www.sasco.co.uk](http://www.sasco.co.uk)



# WARWICK RANGE

## Reusable Medical Polyware

The Warwick Range is a wide range of sterilizable long-life plastic polyware products.

Designed as an alternative to stainless steel products, they are ideal for use in Decontamination Services, Operating Rooms, Clinics, Surgeries and other surgical or healthcare services.

Products are designed for repeated washing, disinfection and sterilization. The Warwick Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>). Products comply with 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

## Warwick Benefits

- Low cost.
  - Robust and hard wearing - guaranteed for 1000 cycles.
  - Easy to clean - blood protein is easily removed from the surface, leaving no residue.
  - 50% lighter than stainless steel, assisting users with storage, lifting and handling.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).
  - Can be used as standalone items or as a support component in surgical sets.
  - Widely used in Sterile Services for packing instrument sets and procedure kits for the autoclave.
- Disinfection sets can be used for manual sterilization in Sterile Services, Clinics, etc.
  - Used non-sterile in hospital wards and clinics.
  - Smooth corners and surfaces, reducing risk of paper rips and the need for re-processing.
  - Autoclave sterilization temperature is 134.6°C + 2°C (generally referred to as 135°C).

## Compartment Trays, Disinfectant and Soaking Containers



Compartment Tray - 3 compartments

Code: CIT1814

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 38

R W S E HP LF



Compartment Tray - 3 compartments

Code: CIT2013

Dimensions: 205x130x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 120

R W S E HP LF



Compartment Tray - 4 compartments

CIT2718

Dimensions: 270x180x41mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Compartment Tray - 4 compartments

Code: LCT2215

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160

R W S E HP LF



Disinfection Set - Soaking container, strainer and hinged lid  
Code: DSSW1000

Dimensions: 260x130x65mm  
Volume: 1litre  
Material: Polypropylene  
Colour: White  
Pack Size: 10

R W S E HP LF



Disinfection Set - Soaking container, strainer and lid  
Code: DSS2711

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer and lid  
Code: DSSW2711

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: White  
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer and vented lid  
Code: DSSV2711

Dimensions: 257x180x128mm  
Volume: 3litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 5

R W S E HP LF



Disinfection Set - Soaking container, strainer, drip tray and lid  
Code: DSS4215

Dimensions: 426x156x138mm  
Volume: 4litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF



Disinfection Set - Soaking container, strainer, drip tray and lid  
Code: DSS8515

Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF



Disinfection Set - 2 x soaking containers, 1 x strainer, 2 x lids and 1 x drip tray  
Code: DSS8515TP


Dimensions: 890x200x145mm  
Volume: 12litre  
Material: Polypropylene  
Colour: Blue  
Pack Size: 1

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free



DIN Style Baskets and Instrument Trays




Instrument Tray - Solid base

**Code: IT1813**

Dimensions: 181x134x54mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 40

R W S E HP LF




Instrument Tray - Solid base

**Code: IT1913**

Dimensions: 195x135x23mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 115

R W S E HP LF




Instrument Tray - Solid base

**Code: T2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Instrument Tray - Solid ribbed base

**Code: IT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Instrument Tray - Solid base

**Code: T3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Instrument Tray - Solid base

**Code: IT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF


**How to order**

To place an order, request a quotation or for general enquiries:

CALL us on  
**+44 (0)1926 422 427**

FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)




Instrument Tray - Solid ribbed base

**Code: IT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Instrument Tray - Solid ribbed base


**Code: IT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12

R W S E HP LF

R Reusable      W Washer Decontamination      S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization      HP Suitable for Hydrogen Peroxide Sterilization      LF Latex Free

DIN Style Baskets and Instrument Trays




Instrument Tray - Mesh base and sides

**Code: MIT1515**

Dimensions: 200x180x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 50

R W S E HP LF




Instrument Tray - Mesh base

**Code: MIT2015**

Dimensions: 200x150x51mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Instrument Tray - Mesh base and sides

**Code: MIT2709**

Dimensions: 275x110x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Instrument Tray - Mesh base

**Code: MIT2715**

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Instrument Tray - Mesh base and sides

**Code: MIT2817**

Dimensions: 280x170x55mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF




Instrument Tray - Mesh base and sides

**Code: MIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 24

R W S E HP LF




Instrument Tray - Mesh base

**Code: MIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12

R W S E HP LF




Instrument Tray - Perforated base

**Code: PIT3025**

Dimensions: 300x250x52mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF




Instrument Tray - Perforated base

**Code: PIT4030**

Dimensions: 424x305x75mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 12

R W S E HP LF




DIN Style Tray - Mesh base and sides

**Code: DINT2525**

Dimensions: 255x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF



DIN Style Tray - Mesh base and sides

**Code: DINT4825**

Dimensions: 480x255x50mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF

R Reusable      W Washer Decontamination      S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization      HP Suitable for Hydrogen Peroxide Sterilization      LF Latex Free

Instrument Tray Lids



Instrument Tray Lid -  
For trays IT1813/CIT1814

Code: ITL181

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -  
For trays IT2015/MIT2015/  
T2015

Code: ITL201

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -  
For trays IT2715/MIT2715

Code: ITL271

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -  
For tray CIT2718

Code: ITL281

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -  
For tray MIT2817

Code: ITL287

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF



Instrument Tray Lid -  
For trays IT3025/MIT3025/  
PIT3025/T3025

Code: ITL302

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 10

R W S E HP LF

How to order

To place an order, request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**  
FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

Gallipots, Lotion and Sponge Bowls



Gallipot

Code: GP40

Dimensions: 40x30mm  
Volume: 35ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 100

R W S E HP LF



Gallipot

Code: GP60

Dimensions: 60x40mm  
Volume: 95ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 300

R W S E HP LF



Gallipot

Code: GP80

Dimensions: 80x60mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 80

R W S E HP LF



Gallipot - Flat base  
Graduated

Code: GPLB60

Dimensions: 60x40mm  
Volume: 50ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 400

R W S E HP LF

Gallipots, Lotion and Sponge Bowls



Gallipot - Flat base  
Graduated

Code: GPLB80

Dimensions: 80x60mm  
Volume: 200ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 200

R W S E HP LF



Lotion Bowl - Flat base  
Graduated

Code: GPLB100

Dimensions: 100x45mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 135

R W S E HP LF



Lotion Bowl - Flat base  
Graduated

Code: GPLB150

Dimensions: 150x70mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Sponge Bowl - Flat base  
Graduated

Code: SB140

Dimensions: 140x65mm  
Volume: 300ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160

R W S E HP LF



Lotion Bowl  
Graduated

Code: LB100

Dimensions: 100x45mm  
Volume: 280ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 200

R W S E HP LF



Lotion Bowl  
Graduated

Code: LB150

Dimensions: 150x70mm  
Volume: 900ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF



Lotion Bowl  
Graduated

Code: LB200

Dimensions: 200x90mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 30

R W S E HP LF



Lotion Bowl  
Graduated

Code: LB250

Dimensions: 250x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF

Jugs



Jug  
Graduated

Code: JG500

Dimensions: N/A  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 50

R W S E HP LF



Jug  
Graduated

Code: JG1000

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 45

R W S E HP LF



Jug  
Graduated


Code: JG2000

Dimensions: N/A  
Volume: 2000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 20

R W S E HP LF



Kidney Dishes, Theatre Bowls and Wash Bowls




Kidney Dish Graduated

**Code: KD150**

Dimensions: 150x35mm  
Volume: 150ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 210

R W S E HP LF




Kidney Dish Graduated

**Code: KD200**

Dimensions: 200x45mm  
Volume: 500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 150

R W S E HP LF




Kidney Dish Graduated

**Code: KD250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 60

R W S E HP LF




Kidney Dish Graduated

**Code: KDY250**

Dimensions: 250x55mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Yellow  
Pack Size: 60

R W S E HP LF




Kidney Dish Graduated

**Code: KD300**

Dimensions: 300x65mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF




Kidney Dish Tray - Shallow kidney shaped tray

**Code: KDT185**

Dimensions: 185x92mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 70

R W S E HP LF




Washbowl Graduated

**Code: WB3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF




Wash Bowl Graduated

**Code: WBN3112**

Dimensions: 310x105mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25

R W S E HP LF




Theatre Bowl Graduated

**Code: WB3513**

Dimensions: 350x120mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Blue  
Pack Size: 25

R W S E HP LF



Wash Bowl Graduated


**Code: WBN3513**

Dimensions: 350x120mm  
Volume: 5000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25

R W S E HP LF

R Reusable      W Washer Decontamination      S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization      HP Suitable for Hydrogen Peroxide Sterilization      LF Latex Free

Medical Boxes, Instrument Boxes & Storage Containers




Instrument Box - Storage box with lid

**Code: IBX2718**

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF




Instrument box - Storage box with vented lid

**Code: IBXV2718**

Dimensions: 257x180x128mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 20

R W S E HP LF



Instrument Box - Storage box with lid

**Code: IBXN290**

Dimensions: 300x224x140mm  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 12

R W S E HP LF


**How to order**

To place an order, request a quotation or for general enquiries:

CALL us on  
**+44 (0)1926 422 427**

FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)




Transportation Box - With hinged lid

**Code: MB2318G**

Dimensions: 230x170x200mm  
Volume: N/A  
Material: ABS  
Colour: Green  
Pack Size: 10

R W S E HP LF



Transportation Box - With hinged lid

**Code: MB2318**

Dimensions: 230x170x200mm  
Volume: N/A  
Material: ABS  
Colour: Orange  
Pack Size: 10

R W S E HP LF

R Reusable      W Washer Decontamination      S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization      HP Suitable for Hydrogen Peroxide Sterilization      LF Latex Free

Scope Containers



Scope Container -  
With silicone dividers and lid

Code: SMT5510

Dimensions: 660x110x66mm  
Volume: N/A  
Material: Polypropylene  
Colour: Lilac  
Pack Size: 1

R W S E HP LF



Silicone Plug -  
For SMT5510

Code: SG5510

Dimensions: N/A  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF

Silicone Protection



Silicone Pad -  
For trays IT/MIT2715 and  
small metal DIN tray  
Code: SP2313

Dimensions: 230x132mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
For trays IT/MIT/PIT/T3025

Code: SP2823

Dimensions: 280x230mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
Fits 1/2 tray IT/MIT/PIT4030

Code: SP2719

Dimensions: 270x190mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
Fits 1/2 tray DINT4825

Code: SP4622

Dimensions: 230x220mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Pad -  
For medium metal DIN tray

Code: SP2324

Dimensions: 235x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF




Silicone Pad -  
For large metal DIN tray

Code: SP4724

Dimensions: 470x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 10

R W S E HP LF



Silicone Roll -  
1 metre length with 10  
divisible sections  
Code: SR10

Dimensions: 100x245mm  
Volume: N/A  
Material: Silicone  
Colour: Pink  
Pack Size: 1

R W S E HP LF

Tray Tags



Tray Tag

Code: TAG-010

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: White  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-015

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Yellow  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-020

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Blue  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-022

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Light Green  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-025

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Orange  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-030

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Red  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-050

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Green  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-055

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Grey  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-060

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Lavender  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-065

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Tan  
Pack Size: 100

R W S E HP LF



Tray Tag

Code: TAG-070

Dimensions: 105x40mm  
Volume: N/A  
Material: Nylon  
Colour: Black  
Pack Size: 100

R W S E HP LF



Stainless Steel Attachment Ring

Code: TR30

Dimensions: 30mm  
Volume: N/A  
Material: Steel  
Colour: Silver  
Pack Size: 100

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free



# STERILITE RANGE

## Single Use Medical Polyware

The SteriLite Range is a wide range of single use plastic polyware products.

Designed to speed up packing and sterilizing procedures. SteriLite is ideal for use in Sterile Pack Manufacturing, Decontamination Services, Operating Rooms and Clinics due to its low cost, product design and performance in hospitals.

The Sterilite Range of products can be processed by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).

Products comply with 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable. Material used is polypropylene unless otherwise stated.

### SteriLite Benefits

- Low cost and disposable.
- No need to invest time and money into retrieval and decontamination after use.
- Clean and ready for sterilization.
- Can be sterilized by steam autoclaving at 135°C, Ethylene Oxide (EO) and Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>).
- Can be used as stand-alone items or as a support (supplementary) component in operating rooms and surgical procedures.
- Widely used in Sterile Services for packing instrument sets and procedure kits.
- Used non-sterile in hospital wards and clinics.



## Compartment Trays and Packing Trays



Packing Tray -  
4 compartments

Code: LCT2215S

Dimensions: 220x150x24mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 160

2 W S E HP LF



Packing Tray

Code: ST4330

Dimensions: 430x330x64mm  
Volume: N/A  
Material: Polypropylene  
Colour: Blue  
Pack Size: 63

2 W S E HP LF

- 2 Single Use
- W Washer Decontamination
- S Suitable for Steam Sterilization
- E Suitable for Ethylene Oxide Sterilization
- HP Suitable for Hydrogen Peroxide Sterilization
- LF Latex Free

Gallipots and Sponge Bowls



			
Gallipot Graduated	Gallipot Graduated	Gallipot Graduated	Gallipot Graduated
Code: GPS60	Code: GPS120	Code: GP100	Code: GP130
Dimensions: 62x30mm	Dimensions: 85x30mm	Dimensions: 100x53mm	Dimensions: 130x66mm
Volume: 60ml	Volume: 120ml	Volume: 250ml	Volume: 500ml
Material: Polypropylene	Material: Polypropylene	Material: Polypropylene	Material: Polypropylene
Colour: Natural	Colour: Natural	Colour: Natural	Colour: Natural
Pack Size: 2000	Pack Size: 1000	Pack Size: 1050	Pack Size: 735
<div>2WSEHP LF</div>	<div>2WSEHP LF</div>	<div>2WSEHP LF</div>	<div>2WSEHP LF</div>



Sponge Bowl Graduated	Sponge Bowl Graduated
Code: SB170	Code: SB200
Dimensions: 170x84mm	Dimensions: 210x90mm
Volume: 1000ml	Volume: 2000ml
Material: Polypropylene	Material: Polypropylene
Colour: Natural	Colour: Natural
Pack Size: 340	Pack Size: 220
<div>2WSEHP LF</div>	<div>2WSEHP LF</div>

- 2 Single Use
- W Washer Decontamination
- S Suitable for Steam Sterilization
- E Suitable for Ethylene Oxide Sterilization
- HP Suitable for Hydrogen Peroxide Sterilization
- LF Latex Free

Jugs and Quivers

	
Jug Graduated	Quiver
Code: J1000	Code: QV30
Dimensions: N/A	Dimensions: 190x65mm
Volume: 1000ml	Volume: N/A
Material: Polypropylene	Material: Polypropylene
Colour: Natural	Colour: Natural
Pack Size: 100	Pack Size: 75
<div>2WSEHP LF</div>	<div>2WSEHP LF</div>

Kidney Dishes and Theatre Bowls

		
Kidney Dish Graduated	Kidney Dish Graduated	Theatre Bowl
Code: KB250	Code: KB300	Code: TB350
Dimensions: 252x100x48mm	Dimensions: 300x140x57mm	Dimensions: 352x116mm
Volume: 750ml	Volume: 1500ml	Volume: 5000ml
Material: Polypropylene	Material: Polypropylene	Material: Polypropylene
Colour: Natural	Colour: Natural	Colour: Blue
Pack Size: 450	Pack Size: 135	Pack Size: 50
<div>2WSEHP LF</div>	<div>2WSEHP LF</div>	<div>2WSEHP LF</div>

How to order

To place an order, request a quotation or for general enquiries:

CALL us on  
+44 (0)1926 422 427  
FAX  
+44 (0)1926 315 405

EMAIL  
enquiries@sasco.co.uk

- 2 Single Use
- W Washer Decontamination
- S Suitable for Steam Sterilization
- E Suitable for Ethylene Oxide Sterilization
- HP Suitable for Hydrogen Peroxide Sterilization
- LF Latex Free



# SPA RANGE

## Nursing & Care, Daily Living Aids and Accessories

The Spa range includes reusable products for personal living, patient care, toileting and bedside comfort.

The Spa Range is designed for use in Hospital Wards, Rehabilitation, Nursing Homes, Community Care and for Personal Use in the Home.

Products comply with BS2588 For Portable Sanitary Products, BS3215 for Reusable Urinals, 93/42/EEC Medical Devices Directive Class 1 and show CE where applicable.

Material used is polypropylene unless otherwise stated.

### Spa Benefits

- Toileting products that preserve patient dignity and independence by empowering patients to toilet themselves with little or no assistance.
- Urinal Bottles specifically designed to the male and female anatomy for a comfortable and safe fit during use.
- Product design reduces risk of cross infection, unpleasant odours and sights.
- Easy to use, clean and disinfect, meaning both patient and carer can be confident that products are safe to re-use.
- Beakers and Drinking Cups are wide, easy to grip, non-scratch, drop-proof and have a nondrip rim.
- Handled beakers are ideal for patients with unsteady or infirm hands.
- Medicine Measures are sturdy, ideal for long-term use and regularly checked by Trading Standards for measurement accuracy and compliance.
- All Spa products can be cleaned in automatic washer disinfectors. Also suitable for hand washing using detergents and disinfectants.



## Bedpans

Bedpan - Adult

Code: BPN

Dimensions: 370x310x90mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 30

R W S E HP LF



Slipper Pan - Bedpan with handle and lid

Code: FPN

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 12

R W S E HP LF

### How to order

To place an order, request a quotation or for general enquiries:

CALL us on  
+44 (0)1926 422 427  
FAX  
+44 (0)1926 315 405

EMAIL  
enquiries@sasco.co.uk



Slipper Pan - Bedpan with handle and lid

Code: FPNW

Dimensions: 449x309x99mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 12

R W S E HP LF



Bedpan - Hospital style bedpan

Code: VP

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 18

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPN

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 15

R W S E HP LF



Bedpan - Hospital bedpan with lid

Code: VPNW

Dimensions: 355x308x108mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 15

R W S E HP LF



Bedpan - European bedpan with handle and lid

Code: DPNF

Dimensions: 420x315x100mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 15

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Bedpans, Commode Pans, Bidets and Bedpan Holders

Slipper Pan -  
Female urinal and  
infants bedpan  
Code: SPN

Dimensions: 335x245mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 25

R W S E HP LF



Commode Pan / Bidet

Code: CMP

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 20

R W S E HP LF



Commode Pan / Bidet

Code: CMPW

Dimensions: 359x338x130mm  
Volume: 3000ml  
Material: Polypropylene  
Colour: White  
Pack Size: 20

R W S E HP LF



Commode Pan -  
Commode pan with lid

Code: CBL

Dimensions: 259x259x144mm  
Volume: 2500ml  
Material: Polypropylene  
Colour: White  
Pack Size: 25

R W S E HP LF



Bedpan Support -  
For pulp disposable midi  
slipper pan 1.3 litres  
Code: PSPS

Dimensions: 470x290x75mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 30

R W S E HP LF



Bedpan Support -  
For pulp disposable bedpan  
2 litres  
Code: PBPS

Dimensions: 370x310x90mm  
Volume: 2000ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 30

R W S E HP LF

How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
**+44 (0)1926 422 427**  
FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

2 Single Use W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Drinking Beakers, Cups, Lids and Denture Cups



Beaker - Feeder/drinking cup  
Graduated

Code: BFC250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Beaker - Feeder/drinking cup  
with handles  
Graduated

Code: BFH250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 140

R W S E HP LF



Lid -  
Narrow Spout for BFC250/  
BFH250

Code: BNS

Dimensions: 4x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Medium Spout for BFC250/  
BFH250

Code: BMS

Dimensions: 7x4mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Wide spout for BFC250  
& BFH250

Code: BWS

Dimensions: 13x8mm spout  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 300

R W S E HP LF



Lid -  
Rim feeder for beakers BFC250  
& BFH250

Code: BRF

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: Natural  
Pack Size: 200

R W S E HP LF



Drinking Cup -  
With handles, spout lid and  
feeder lid

Code: GPMC6250

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: Cream  
Pack Size: 90

R W S E HP LF



Denture Cup -  
With strainer and lid

Code: DNT

Dimensions: 90x80x60mm  
Volume: N/A  
Material: Polypropylene  
Colour: Green  
Pack Size: 75

R W S E HP LF



Beaker Set - Feeder/drinking  
cup with handle, spout lid and  
feeder lid  
Graduated

Code: HDMS

Dimensions: N/A  
Volume: 300ml  
Material: SAN  
Colour: Clear  
Pack Size: 40

R W S E HP LF



Teapot Style Beaker

Code: TFB

Dimensions: N/A  
Volume: 250ml  
Material: Polypropylene  
Colour: White  
Pack Size: 45

R W S E HP LF



Clear Glass Style Beaker

Code: TBC

Dimensions: N/A  
Volume: 280ml  
Material: SAN  
Colour: Clear  
Pack Size: 120

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free



Patient Trays, Kidney Dishes and Vomit Bowls



Dispensing/Dressing Tray

Code: ITC2715

Dimensions: 270x150x30mm  
Volume: N/A  
Material: Polypropylene  
Colour: Green  
Pack Size: 60

R W S E HP LF



Kidney Dish  
Graduated

Code: KDC250

Dimensions: 250x130mm  
Volume: 750ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 60

R W S E HP LF



Vomit Bowl - With handle  
Graduated

Code: VB300

Dimensions: 300x144mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 25

R W S E HP LF

R Reusable W Washer Decontamination S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization HP Suitable for Hydrogen Peroxide Sterilization LF Latex Free

Jugs and Jug Sets



Jug  
Graduated

Code: JG500

Dimensions: N/A  
Volume: 500ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 50

R W S E HP LF



Jug  
Graduated

Code: JG1000

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 45

R W S E HP LF



Jug with lid  
Graduated

Code: JGL1000

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 40

R W S E HP LF



Jug Set - Includes tray, carafe  
with lid and beaker  
Graduated

Code: LST

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: N/A  
Pack Size: 20

R W S E HP LF

How to order

To place an order,  
request a quotation or for  
general enquiries:

CALL us on  
+44 (0)1926 422 427  
FAX  
+44 (0)1926 315 405

EMAIL  
enquiries@sasco.co.uk

Medicine Measures



Medicine Measure -  
Graduated

Code: MM50

Dimensions: N/A  
Volume: 50ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000

R W S E HP LF



Medicine Measure -  
Graduated

Code: MM60

Dimensions: N/A  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000

R W S E HP LF




Medicine Measure -  
Graduated black lines

Code: MM60BL

Dimensions: N/A  
Volume: 60ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 1000

R W S E HP LF

Urinal Bottles and Urinal Pans




Slipper Pan - Female urinal

**Code: SPN**

Dimensions: 335x245mm  
Volume: 1500ml  
Material: Polypropylene  
Colour: Green  
Pack Size: 25

R W S E HP LF




Urinal Bottle - Female urinal with handle

**Code: UF1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 30

R W S E HP LF




Urinal Pan - Female Urinal

**Code: PUF1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: White  
Pack Size: 20

R W S E HP LF




Urinal Bottle - Male urinal with handle Graduated

**Code: UMA1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 40

R W S E HP LF




Urinal Bottle - Male urinal with handle, hospital style Graduated

**Code: UBS1000**

Dimensions: N/A  
Volume: 1200ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 45

R W S E HP LF




Urinal Bottle - Male urinal with handle

**Code: UML1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polyethylene  
Colour: Natural  
Pack Size: 50

R W S E HP LF




Urinal Bottle - Male urinal, hospital style Graduated

**Code: UM1000**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 25

R W S E HP LF




Urinal Bottle - Male urinal, hospital style Graduated

**Code: UM650**

Dimensions: N/A  
Volume: 650ml  
Material: Polypropylene  
Colour: Natural  
Pack Size: 30

R W S E HP LF




Non Leak Valve - For use with UMA1000 & UML1000

**Code: UMT**

Dimensions: N/A  
Volume: N/A  
Material: Polypropylene  
Colour: White  
Pack Size: 25

R W S E HP LF




Urinal Bottle with UMT -Male urinal bottle with non leak valve Graduated

**Code: USP1200**

Dimensions: N/A  
Volume: 1000ml  
Material: Polypropylene and Polyethylene  
Colour: N/A  
Pack Size: 10

R W S E HP LF



Urinal Bottle Holder - For UBS1000, UML1000 and UM1000

**Code: UBH2002**

Dimensions: N/A  
Volume: 1200ml  
Material: Plastic Coated Wire  
Colour: White  
Pack Size: 10

R W S E HP LF

**How to order**

To place an order, request a quotation or for general enquiries:

CALL us on  
**+44 (0)1926 422 427**  
FAX  
**+44 (0)1926 315 405**

EMAIL  
[enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

R Reusable      W Washer Decontamination      S Suitable for Steam Sterilization  
E Suitable for Ethylene Oxide Sterilization      HP Suitable for Hydrogen Peroxide Sterilization      LF Latex Free





Warwick House, Heathcote Way, Warwick, CV34 6TE  
T: +00 44 (0) 1926 422 427  
F: +00 44 (0) 1926 315 405  
E: [enquiries@sasco.co.uk](mailto:enquiries@sasco.co.uk)

**[www.sasco.co.uk](http://www.sasco.co.uk)**